

# VISUALIZATION PLATFORM


## Sensors Data Analysis For Assisting Elderly People

### ACCESS

- 1 Through the website
- 2 First login with your username and password

#### S5 - Visualization Platform

Sensors data analysis platform for assisting elderly people



### OBJECTS & DATA IMPORTS

#### Objects

This visualization platform aims to give you, as doctor of patients, the results and analysis of the sensors placed in you patient smart home.

Actually the website has three main visualization pages :

- ✓ **ADLs** : giving you information of the **performed activities of daily living**
- ✓ **Dependency** : giving you details about the **correlations** between the different performed activities
- ✓ **Sensors** : giving you information about the **geographic positions** and the **raw sensors data**

#### Data imports


##### ➤ Activities Dataset

On the navigation bar, click on “Import Dataset” as follow :

An import system page opens and enables you to upload the activities dataset of you patient.

Be careful to respect the format of the file that as to be according to the following scheme “**date time activity status**” as shown below :


```
2008-02-25 00:19:32 brush teeth start
2008-02-25 00:21:23 brush teeth end
2008-02-25 00:22:46 go to bed start
```



##### ➤ Home and Sensors Dataset

On the navigation bar, click on “Import HomeDataset” as follow :

An import system page opens and enables you to upload **three files** : the first being the map of you patient home (**image file**), the second informing about the sensors coordinates (in pixel) With the format “**sensorName[tab]width[tab]height**”, the last one giving the sensors data (with header) with format “**startDatetime[tab]endDatetime[end]SensorName**”



| 3 | start time          | End time            | Sensor            |
|---|---------------------|---------------------|-------------------|
|   | 2008-02-25 00:20:14 | 2008-02-25 00:22:57 | Hall-Bedroom door |

| 2 | frontdoor   | 519 | 716 |
|---|-------------|-----|-----|
|   | Dishwasher  | 710 | 130 |
|   | ToiletFlush | 253 | 562 |



### SELECTIONS

➤ Patient Selection

➤ Dates Selection

PatientName

From: 25/Feb/2008

To: 25/Feb/2008

Possible from the navigation bar

# Visualization Platform : ADLs page (default page)

## First Graph : Status Table



### Example of Status Table :

| Event       | LastTime  |
|-------------|---|
| Go to bed   |  |
| Use toilet  |  |
| Brush teeth |  |

### Status Table

This table shows the status for each activity.



### Measures

- 1 **Last time green dot** : the last time the activity has been performed is less than 0.5 days.
- 2 **Last time orange dot** : the last time the activity has been performed is between 0.5 and one day.
- 3 **Last time red dot** : the last time the activity has been performed is more than one day.

### Filters

- 4 **Patient** : selection of the patient from the navigation bar. Result given for this patient.
- 5 **Dates** : selection of a date of start and a date of end to consider activities data.  
Not useful here as it is compared with the current date time.



### How to read the graph?

This graph provides information about the last time an activity has been performed by the patient. If the dot is red, it may mean that an anomaly has occurred.

# Visualization Platform : ADLs page (default page)


## Second Graph : Time Consuming Activities Pie Chart



### Example of a pie chart and its details table:

#### Time consuming activities

This pie chart shows the time consuming activities of the patient. For more details, please refer to the table below.



#### Details

This table details the above pie chart. For each event, you can find the frequency of occurrences and also the proportion of time it represents compared to the period you consider (whether day, week or month).

| Event             | Frequency | Duration |
|-------------------|-----------|----------|
| Brush teeth       | 3         | 0.06 %   |
| Get drink         | 3         | 0.04 %   |
| Go to bed         | 3         | 19.15 %  |
| Leave house       | 2         | 13.77 %  |
| Prepare breakfast | 2         | 0.14 %   |
| Prepare dinner    | 1         | 1.48 %   |
| Take shower       | 2         | 0.36 %   |
| Use toilet        | 11        | 0.03 %   |
| Unknown           | 1         | 64.98 %  |



### Measures

- 1 **Event** : the performed activity during the considered period
- 2 **Frequency** : the occurrences of the performed activity during the given period.
- 3 **Duration** : the proportion of time it represents regarding to the selected period (in percents).

### Filters

- 4 **Patient** : selection of the patient from the navigation bar. Result given for this patient.
- 5 **Dates** : selection of a date of start and a date of end to consider activities data.




### How to read the graph?

These graph shows you how the time is spent during a selected period. The table completes the pie chart by given all the results numbers included the smallest ones. You will then be able to check how often an activity is realized by the patient in a given time.



Example of a circular heat chart:



Measures

- 1 **Hour of the day** : from midnight to 11pm.
- 2 **Activity**: the performed activity during the considered period
- 3 **Frequency Gradient Color** : the occurrences of the performed activity at the given hour of the day during the considered period


Filters

- 4 **Patient** : selection of the patient from the navigation bar. Result given for this patient.
- 5 **Dates** : selection of a date of start and a date of end to consider activities data.

How to read the graph?

The graph enables you to check the habits of your patient and deduce a potential typical day. Indeed you can for instance see when the patient is used to going to bed or eating.

Example of a chord diagram:



Measures


- 1 **Relation between activities branch** : activities directly related (one happens just after the other) and frequency of this relation thanks to the length of the branch
- 2 **Activity**: the performed activity during the considered period

Filters

- 3 **Patient** : selection of the patient from the navigation bar. Result given for this patient.
- 4 **Dates** : selection of a date of start and a date of end to consider activities data.

How to read the graph?

First you can focus on the relation of one activity by clicking On it. It will fade the other relations away. Then the length of the branch on the considered activity side represents how often it is realized after the linked activity and, the same way, the length of the same branch on the linked activity side displays the times the considered activity is followed by this one (inversed relation)




# Visualization Platform : Dependency page

## Second Graph : Gantt Chart

GRAPH : GANTT CHART



### Example of a gantt chart:



PARAMETERS



### Measures

- 1 **Time** : time of day and thus duration during which activity is performed.
- 2 **Activity**: the performed activity during the considered period

### Filters

- 3 **Patient** : selection of the patient from the navigation bar. Result given for this patient.
- 4 **Dates** : selection of a date of start and a date of end to consider activities data.

UTILIZATION




### How to read the graph?

This Gantt graph enables you to see the sequences of the activities during a given period. Each colored segments represents the duration of the activity. We can thus read the date and time of start and end of this activity too.



Example of a tree map:



Tree Map

This tree map displays the frequency of occurrences of an event following another one.



Measures

- 1 **Activity**: the performed activity during the considered period given by a specific color.
- 2 **Next Activity Frequency**: the occurrences of the activity following the considered activity given by the size of the rectangle.

Filters

- 3 **Patient** : selection of the patient from the navigation bar. Result given for this patient.
- 4 **Dates** : selection of a date of start and a date of end to consider activities data.




How to read the graph?

Each color represents a performed activity during the considered period. The rectangles represents the sequence frequency of another activity following this one. The sum of all the rectangles of the same color give logically the frequency of the corresponding activity.

On a click or hover in a rectangle we can see the frequency value, its proportion (in percent) regarding the parent activity and finally its proportion regarding all rectangles.



Example of a heat map and its details table:



Heat Map

This map gives the frequency of the geographic points of the patients in a given time.

Details

This table details the above heat map. For each sensor, you can find the frequency of occurrences and the average duration of the sensor activation.

| Sensor             | Frequency | AvgDuration |
|--------------------|-----------|-------------|
| Hall-Toilet door   | 3         | 96 sec.     |
| Hall-Bathroom door | 7         | 5 sec.      |
| Plates cupboard    | 2         | 4 sec.      |
| Frontdoor          | 1         | 5 sec.      |
| ToiletFlush        | 3         | 1 sec.      |
| Freezer            | 3         | 9 sec.      |
| Washingmachine     | 3         | 39 sec.     |
| Hall-Bedroom door  | 4         | 1016 sec.   |



Measures

- 1 **Sensor** : the activated sensor during the considered period.
- 2 **Frequency** : the occurrences of the sensor activation given by a gradient color on its geographic point on the map.
- 3 **AvgDuration** : the mean time the sensor is activated during the considered period.

Filters

- 4 **Patient** : selection of the patient from the navigation bar. Result given for this patient.
- 5 **Dates** : selection of a date of start and a date of end to consider activities data.



How to read the graph?

The heat map will show you the most frequent position taken by the patient at a given period. The frequency is represented by a gradient color going from the blue (cold color) for few passages to yellow for medium passages to red for very frequent positions. A details table completes the heat map graph by giving for each sensor the frequency exact value and also providing the average time the sensor stays activated during the considered period.